



## 5 oz. Snap Cap with Lid

### Product Code:

0450-1110 (Non-Sterile)

0450-1111 (Sterile)

Product is individually wrapped



Filled with dye to illustrate details

### General Description

Color: Transparent

Body Diameter: 73.66 mm

Body Height: 66.29 mm

- Snap-on lid creates secure seal
- Shatter-resistant, Transparent body allows for easy viewing
- Graduated in .5oz and 10mL increments
- Can also be Sterilized

### Material:

Polypropylene, FDA Certified

Country of Origin: China

Product is recyclable

### Packing Details

Case: 100 pc/case




Packaged: Individually

Case Dimensions: 15.25 x 15 x 14.50"

Case Weight: 6.1 lbs

Cube Dimension: 1.92

Cases per Skid: 30



Manufactured by Oak Ridge Products, LLC  
FDA Manufacturing Facility Registration  
#1419823 & 3006060667

OAK RIDGE PRODUCTS, LLC.  
ALL RIGHTS RESERVED

Dimensions listed are in millimeters

US: 888-650-7444 INT'L: 001-815-363-4700  
WWW.OAKRIDGEPRODUCTS.COM REV. 1/10/14