



8 oz. Snap Cap with Lid

Product Code:

0480-1110 (Non-Sterile)

0480-1111 (Sterile)

Product is individually wrapped



Filled with dye to illustrate details

General Description

Color: Transparent

Body Diameter: 73.66 mm

Body Height: 100.51 mm

• Snap-on lid creates secure seal

• Shatter-resistant, Transparent body allows for easy viewing

• Graduated in .5oz and 10mL increments

• Can also be Sterilized

Material, Packing

Polypropylene, FDA Certified

Country of Origin: China

Product is recyclable

Packing Details

Case: 100 pc/case



Packaged: Individually

Case Dimensions: 15.25 x 15 x 19.25"

Case Weight: 7.4 lbs

Cube Dimension: 2.55

Cases per Skid: 30



Manufactured by Oak Ridge Products, LLC

FDA Manufacturing Facility Registration

#1419823 & 3006060667

Dimensions listed are in millimeters

OAK RIDGE PRODUCTS, LLC.
ALL RIGHTS RESERVED

US: 888-650-7444 INT'L: 001-815-363-4700
WWW.OAKRIDGEPRODUCTS.COM

REV.1/10/14