



5 oz. Snap Cap without Lid

Product Code:

0450-1100 (Non-Sterile)

0450-1101 (Sterile)

General Description

Color: Transparent

Body Diameter: 73.66 mm

Body Height: 66.29 mm

- Snap-on lid creates secure seal
- Shatter-resistant, Transparent body allows for easy viewing
- Graduated in .5oz and 10mL increments
- Lid sold separately; Product # 0501-3100



Filled with dye to illustrate details

Material:

Polypropylene, FDA Certified

Country of Origin: China

Product is recyclable

Packing Details



Case: 500 pc/case

Packaged: 20 Bagged

Case Dimensions: 14.50 x 14.50 x 12.75"

Case Weight: 14.1 lbs

Cube Dimension: 1.55



Manufactured by Oak Ridge Products, LLC

FDA Manufacturing Facility Registration

#1419823 & 3006060667

Dimensions listed are in millimeters