

## 1oz Medicine Cups

Product Code: MC79



## **General Description**

## **Graduations:**

ounces: 1/8, 1/4, 1/2, 3/4, 1

drams: 1 - 8

tablespoons: ½, 1, 2,

teaspoons: ½, 1

dssp: 1,2

cc: 2 ½, 5, 7 ½, 10, 15, 20, 25, 30


milliliters: 2 ½, 5, 7 ½, 10, 15, 20, 25, 30


## Material:

Polypropylene, FDA Certified Country of Origin: China Product is recyclable

Case: 5000 pc/case
Packaged: 50 sleeves of 100
Case Dimensions: 18 x 9 x 12.5
Case Weight: 22.35 lbs

Cube Dimension: 1.19 Cases per Skid: 50





Manufactured by Oak Ridge Products, LLC FDA Manufacturing Facility Registration #1419823 & 3006060667

Product dimensions listed are in millimeters

OAK RIDGE PRODUCTS, LLC. ALL RIGHTS RESERVED

US:888-650-7444 INT'L:001-815-363-4700 WWW.OAKRIDGEPRODUCTS.COM REV. 1/10/14