



8 oz. Snap Cap without Lid

Product Code:

0480-1100 (Non-Sterile)

0480-1101 (Sterile)



Filled with dye to illustrate details

General Description

Color: Transparent

Body Diameter: 73.66 mm

Body Height: 100.51 mm

- Snap-on lid creates secure seal
- Shatter-resistant, Transparent body allows for easy viewing
- Graduated in .5 oz and 10 mL increments
- Lid sold separately; Product # 0501-3100

Material:

Polypropylene, FDA Certified

Country of Origin: China

Product is recyclable

Packing Details

Case: 500 pc/case



Packaged: 20 Bagged

Case Dimensions: 14.50 x 14.50 x 14"

Case Weight: 20.6 lbs

Cube Dimension: 1.70

Cases per Skid: 30



Manufactured by Oak Ridge Products, LLC

FDA Manufacturing Facility Registration

#1419823 & 3006060667

Dimensions listed are in millimeters

OAK RIDGE PRODUCTS, LLC.
ALL RIGHTS RESERVED

US: 888-650-7444 INT'L: 001-815-363-4700
WWW.OAKRIDGEPRODUCTS.COM

REV. 1/10/14