



24 Hour Collector

Product Code:

0803-2100 (Natural)

0803-2700 (Amber)

General Description

- Graduated in 100 mL increments up to 3000 mL

Color: Natural and Amber

Body Diameter: 115.80 mm

Body Height: 291.02 mm

Material:

HDPE, FDA Certified

Country of Origin: USA

Product is recyclable



Case: 40 pc/case


Packaged: Individual

Case Dimensions: 23.50 x 19.25 x 23.75"

Case Weight: 13.2 lbs.

Cube Dimension: 6.22

Cases per Skid: 12



Manufactured by Oak Ridge Products, LLC
FDA Manufacturing Facility Registration
#1419823 & 3006060667

Dimensions listed are in millimeters