



5 ml Transport Tube

Product Code: 0205-0100



Filled with dye to illustrate details

General Description



- Round bottom tube w/ screw tight cap
- For general sample collection, storage, transport, & processing
- Leak proof depending on closure system and usage

Material:

Polypropylene, FDA Certified
 Country of Origin: China
Product is recyclable

Packing Details

Case:	1000 pc/ctn
Packaged:	2 bags of 500
Case Dimensions:	10.50 x 11 x 8.75"
Case Weight:	9.3 lbs
Cube Dimension:	0.58



Manufactured by Oak Ridge Products, LLC
 FDA Manufacturing Facility Registration
 #1419823 & 3006060667

Product dimensions listed are in millimeters

OAK RIDGE PRODUCTS, LLC.
 ALL RIGHTS RESERVED

US:888-650-7444 INT'L:001-815-363-4700
 WWW.OAKRIDGEPRODUCTS.COM

REV. 1/10/14