



25 ml Vial w/ Colored Cap

Product Code:

0625-0300 (White cap)

0625-0400 (Blue cap)



General Description

- Ideal for general purpose storage and transportation
- Conical bottom for centrifuge
- Cap & vial packaged separately

Material:


Polypropylene Body, FDA Certified

HDPE Cap, FDA Certified

Country of Origin: China

Product is recyclable

Case:	500 pc/case
Packaged:	500 Bagged
Case Dimensions:	12.75 x 13 x 18
Case Weight:	11.2 lbs.
Cube Dimension:	1.63
Cases per Skid:	36



Manufactured by Oak Ridge Products, LLC

FDA Manufacturing Facility Registration
#1419823 & 3006060667

Product dimensions listed are in millimeters

OAK RIDGE PRODUCTS, LLC.
ALL RIGHTS RESERVED

US:888-650-7444 INT'L:001-815-363-4700
WWW.OAKRIDGEPRODUCTS.COM

REV. 1/10/14