



3 Gallon R Style Sharps Container

Product Code: 0330-150R

General Description

- Designed for Multi purpose use.
- Rotary Lid maximizes usable container volume for larger sharps
- Shatter & Puncture-resistant
- Locks for final disposal

Material:


High Impact Polypropylene

Country of Origin: China

Product is recyclable

Packing Details

Case:	10 each/case
Packaged:	Bulk, Nested
Case Dimensions:	14" x 12" x 22"
Case Weight:	16lbs
Cube Dimension:	2.14
Cases per Skid:	12



Manufactured by Oak Ridge Products, LLC
FDA Manufacturing Facility Registration
#1419823 & 3006060667

Product dimensions listed are in inches

OAK RIDGE PRODUCTS, LLC.
ALL RIGHTS RESERVED

US:888-650-7444 INT'L:001-815-363-4700
WWW.OAKRIDGEPRODUCTS.COM

REV. 1/10/14