



2 Gallon Sharps Container

Product Code: 0320-1500

General Description

- Designed for use in OR, Surgery areas. For use by preofessional health caregivers
- Puncture-resistant, transparent body allows for easy viewing

Material:




High Impact Polypropylene

Country of Origin: China

Product is recyclable

Packing Details

Case:	20 each/case
Packaged:	Each
Case Dimensions:	18.5 x 14.125 x 18.5"
Case Weight:	22.6 lbs
Cube Dimension:	2.51
Cases per Skid:	30



Manufactured by Oak Ridge Products, LLC
 FDA Manufacturing Facility Registration
 #1419823 & 3006060667

Product dimensions listed are in millimeters

OAK RIDGE PRODUCTS, LLC.
 ALL RIGHTS RESERVED

US:888-650-7444 INT'L:001-815-363-4700
 WWW.OAKRIDGEPRODUCTS.COM

REV. 1/10/14