

## 2 Gallon Flip Lid Sharps

Product Code: 0320-150F



## **General Description**

- •Designed for Multi purpose use.
- •Flip Lid maximizes usable container volume for larger sharps
- Shatter & Puncture-resistant
- Locks for final disposal

## Material:

High Impact Polypropylene Country of Origin: China


FDA Approved


## **Packing Details**

Case: 20 each/case Packaged: Bulk, Nested

Case Dimensions: 21 x 20.625 x 14.625"

Case Weight: 21 lbs
Cube Dimension: 3.67
Cases per Skid: 20





Manufactured by Oak Ridge Products, LLC FDA Manufacturing Facility Registration #1419823 & 3006060667

Product dimensions listed are in inches